



# KCNIP4 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA750882LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6PIL6
<b>Immunogen</b>	Recombinant Human Kv channel-interacting protein 4 protein (1-250AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Kv channel-interacting protein 4 (KCHIP4) (A-type potassium channel modulatory protein 4) (Calsenilin-like protein) (Potassium channel-interacting protein 4), KCNIP4, CALP KCHIP4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	KCNIP4