



GFRAL Antibody, Biotin conjugated

Product Code	CSB-PA751020LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXV0
Immunogen	Recombinant Human GDNF family receptor alpha-like protein (19-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDNF family receptor alpha-like, GFRAL, C6orf144
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GFRAL