





FBXW12 Antibody

Product Code	CSB-PA751062LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6X9E4
Immunogen	Recombinant Human F-box/WD repeat-containing protein 12 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 12 (F-box and WD-40 domain-containing protein 12) (F-box only protein 35), FBXW12, FBW12 FBXO12 FBXO35
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW12

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA751062LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA751062LA01HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.