



FBXW12 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA751062LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6X9E4 |
| Immunogen | Recombinant Human F-box/WD repeat-containing protein 12 protein (1-250AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | F-box/WD repeat-containing protein 12 (F-box and WD-40 domain-containing protein 12) (F-box only protein 35), FBXW12, FBW12 FBXO12 FBXO35 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | FBXW12 |