



GBP6 Antibody, HRP conjugated

Product Code	CSB-PA751154LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZN66
Immunogen	Recombinant Human Guanylate-binding protein 6 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylate-binding protein 6 (GTP-binding protein 6) (GBP-6) (Guanine nucleotide-binding protein 6), GBP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GBP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.