



DPH6 Antibody, Biotin conjugated

Product Code	CSB-PA752102LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8W6
Immunogen	Recombinant Human Diphthine--ammonia ligase protein (8-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine--ammonia ligase (EC 6.3.1.14) (ATP-binding domain-containing protein 4) (Diphthamide synthase) (Diphthamide synthetase) (Protein DPH6 homolog), DPH6, ATPBD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DPH6