



DPH6 Antibody, Biotin conjugated

Product Code	CSB-PA752102LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8W6
Immunogen	Recombinant Human Diphthine--ammonia ligase protein (8-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine--ammonia ligase (EC 6.3.1.14) (ATP-binding domain-containing protein 4) (Diphthamide synthase) (Diphthamide synthetase) (Protein DPH6 homolog), DPH6, ATPBD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DPH6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.