



clfB Antibody, Biotin conjugated

Product Code	CSB-PA753442HD01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GDH2
Immunogen	Recombinant Staphylococcus aureus Clumping factor B protein (45-542AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Clumping factor B (Fibrinogen receptor B) (Fibrinogen-binding protein B), clfB
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Research Area	Others
Target Names	clfB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.