



FAM109B Antibody

Product Code	CSB-PA753688LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ICB4
Immunogen	Recombinant Human Sesquipedalian-2 protein (1-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sesquipedalian-2 (Ses2) (27 kDa inositol polyphosphate phosphatase interacting protein B) (IPIP27B), FAM109B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM109B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.