



FAM109B Antibody, Biotin conjugated

Product Code	CSB-PA753688LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ICB4
Immunogen	Recombinant Human Sesquipedalian-2 protein (1-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sesquipedalian-2 (Ses2) (27 kDa inositol polyphosphate phosphatase interacting protein B) (IPIP27B), FAM109B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM109B