



ETFRF1 Antibody, FITC conjugated

Product Code	CSB-PA753708LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IPR1
Immunogen	Recombinant Human Electron transfer flavoprotein regulatory factor 1 protein (1-88AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Electron transfer flavoprotein regulatory factor 1 (LYR motif-containing protein 5), ETFRF1, LYRM5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ETFRF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.