



FAM166A Antibody, Biotin conjugated

Product Code	CSB-PA753727LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6J272
Immunogen	Recombinant Human Protein FAM166A protein (1-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM166A, FAM166A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM166A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.