




NANOGP8 Antibody

Product Code	CSB-PA754125LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSW7
Immunogen	Recombinant Human Homeobox protein NANOGP8 protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein NANOGP8, NANOGP8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NANOGP8

Image



Western Blot

Positive WB detected in: Rat liver tissue, Mouse liver tissue

All lanes: NANOGP8 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.