



NANOGP8 Antibody, HRP conjugated

Product Code	CSB-PA754125LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSW7
Immunogen	Recombinant Human Homeobox protein NANOGP8 protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein NANOGP8, NANOGP8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NANOGP8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.