



REEP3 Antibody, Biotin conjugated

Product Code	CSB-PA754137LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUK4
Immunogen	Recombinant Human Receptor expression-enhancing protein 3 protein (80-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor expression-enhancing protein 3, REEP3, C10orf74
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	REEP3