





FAM117B Antibody

Product Code	CSB-PA754170LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L5
Immunogen	Recombinant Human Protein FAM117B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM117B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein), FAM117B, ALS2CR13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM117B

Image



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA754170LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA754170LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.