



FAM117B Antibody, HRP conjugated

Product Code	CSB-PA754170LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L5
Immunogen	Recombinant Human Protein FAM117B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM117B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein), FAM117B, ALS2CR13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM117B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.