



FAM117B Antibody, HRP conjugated

Product Code	CSB-PA754170LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L5
Immunogen	Recombinant Human Protein FAM117B protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM117B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein), FAM117B, ALS2CR13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM117B