



DCUN1D2 Antibody

Product Code	CSB-PA754270LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PH85
Immunogen	Recombinant Human DCN1-like protein 2 protein (1-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DCN1-like protein 2 (DCUN1 domain-containing protein 2) (Defective in cullin neddylation protein 1-like protein 2), DCUN1D2, C13orf17 DCUN1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DCUN1D2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.