



FAM171B Antibody, HRP conjugated

Product Code	CSB-PA757601LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P995
Immunogen	Recombinant Human Protein FAM171B protein (59-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM171B, FAM171B, KIAA1946
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FAM171B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.