



GJB7 Antibody, Biotin conjugated

Product Code	CSB-PA757637LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PEY0
Immunogen	Recombinant Human Gap junction beta-7 protein (41-75AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gap junction beta-7 protein (Connexin-25) (Cx25), GJB7, CX25
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GJB7