



# KCTD18 Antibody

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA757654LA01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q6PI47   |
| Immunogen           | Recombinant Human BTB/POZ domain-containing protein KCTD18 protein (1-200AA)     |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Non-conjugated   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | BTB/POZ domain-containing protein KCTD18, KCTD18                                 |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Signal Transduction  |
| Target Names        | KCTD18   |