



KCTD18 Antibody

Product Code	CSB-PA757654LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PI47
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD18 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing protein KCTD18, KCTD18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCTD18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.