



TPTE2 Antibody, HRP conjugated

Product Code	CSB-PA757850LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XPS3
Immunogen	Recombinant Human Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase TPTE2 protein (167-411AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase TPTE2 (EC 3.1.3.67) (Lipid phosphatase TPIP) (TPTE and PTEN homologous inositol lipid phosphatase), TPTE2, TPIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TPTE2