




inIA Antibody

Product Code	CSB-PA758331HA01AAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q723K6
Immunogen	Recombinant <i>Listeria monocytogenes</i> serotype 4b Internalin-A protein (562-770AA)
Raised In	Rabbit
Species Reactivity	<i>Listeria monocytogenes</i> serotype 4b
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Internalin-A, inIA
Immunogen Species	<i>Listeria monocytogenes</i> serotype 4b
Research Area	Others
Target Names	inIA

Image



Western blot
 All lanes: inIA antibody at 10µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 87 kDa
 Observed band size: 87 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.