



inIA Antibody, FITC conjugated

Product Code	CSB-PA758331HC01AAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q723K6
Immunogen	Recombinant Listeria monocytogenes serotype 4b Internalin-A protein (562-770AA)
Raised In	Rabbit
Species Reactivity	Listeria monocytogenes serotype 4b
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Internalin-A, inIA
Immunogen Species	Listeria monocytogenes serotype 4b (strain F2365)
Research Area	Others
Target Names	inIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.