





ASTL Antibody

Product Code	CSB-PA760708LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6HA08
Immunogen	Recombinant Human Astacin-like metalloendopeptidase protein (38-337AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Astacin-like metalloendopeptidase (EC 3.4.-.-) (Oocyte astacin) (Ovastacin), ASTL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ASTL

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA760708LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA760708LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.