



KRT80 Antibody, Biotin conjugated

Product Code	CSB-PA760849LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6KB66
Immunogen	Recombinant Human Keratin, type II cytoskeletal 80 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 80 (Cytokeratin-80) (CK-80) (Keratin-80) (K80) (Type-II keratin Kb20), KRT80, KB20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT80
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.