



BCOR Antibody

Product Code	CSB-PA761519LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6W2J9
Immunogen	Recombinant Human BCL-6 corepressor protein (1569-1755AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BCL-6 corepressor (BCoR), BCOR, KIAA1575
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	BCOR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.