



HACD2 Antibody, HRP conjugated

Product Code	CSB-PA761567LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Y1H2
Immunogen	Recombinant Human Very-long-chain (3R)-3-hydroxyacyl-CoA dehydratase 2 protein (7-41AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Very-long-chain (3R)-3-hydroxyacyl-CoA dehydratase 2 (EC 4.2.1.134) (3-hydroxyacyl-CoA dehydratase 2) (HACD2) (Protein-tyrosine phosphatase-like member B), HACD2, PTPLB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HACD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.