



FGD6 Antibody, Biotin conjugated

Product Code	CSB-PA761649ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZV73
Immunogen	Recombinant Human FYVE, RhoGEF and PH domain-containing protein 6 protein (1138-1352AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FYVE, RhoGEF and PH domain-containing protein 6 (Zinc finger FYVE domain-containing protein 24), FGD6, KIAA1362 ZFYVE24
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FGD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.