





PAXIP1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA761652LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6ZW49 |
| Immunogen | Recombinant Human PAX-interacting protein 1 protein (613-812AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | PAX-interacting protein 1 (PAX transactivation activation domain-interacting protein), PAXIP1, PAXIP1L PTIP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PAXIP1 |

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA761652LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA761652LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.