



NAPEPLD Antibody, Biotin conjugated

Product Code	CSB-PA764247LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IQ20
Immunogen	Recombinant Human N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D (N-acyl phosphatidylethanolamine phospholipase D) (NAPE-PLD) (NAPE-hydrolyzing phospholipase D) (EC 3.1.4.54), NAPEPLD, C7orf18
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NAPEPLD