



MYORG Antibody, FITC conjugated

Product Code	CSB-PA764657LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSJ0
Immunogen	Recombinant Human Myogenesis-regulating glycosidase protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myogenesis-regulating glycosidase (EC 3.2.1.-) (Uncharacterized family 31 glucosidase KIAA1161), MYORG, KIAA1161
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MYORG