



MYORG Antibody, Biotin conjugated

Product Code	CSB-PA764657LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSJ0
Immunogen	Recombinant Human Myogenesis-regulating glycosidase protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myogenesis-regulating glycosidase (EC 3.2.1.-) (Uncharacterized family 31 glucosidase KIAA1161), MYORG, KIAA1161
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MYORG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.