



MED27 Antibody, Biotin conjugated

Product Code	CSB-PA764716LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P2C8
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 27 protein (39-108AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 27 (Cofactor required for Sp1 transcriptional activation subunit 8) (CRSP complex subunit 8) (Mediator complex subunit 27) (P37 TRAP/SMCC/PC2 subunit) (Transcriptional coactivator CRSP34), MED27, CRSP34 CRSP8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MED27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.