



PGM2L1 Antibody, HRP conjugated

Product Code	CSB-PA764776LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PCE3
Immunogen	Recombinant Human Glucose 1,6-bisphosphate synthase protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose 1,6-bisphosphate synthase (EC 2.7.1.106) (PMMLP) (Phosphoglucomutase-2-like 1), PGM2L1, BM32A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PGM2L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.