



TRAT1 Antibody, Biotin conjugated

Product Code	CSB-PA764808LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PIZ9
Immunogen	Recombinant Human T-cell receptor-associated transmembrane adapter 1 protein (29-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30), TRAT1, TCRIM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TRAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.