



SPATA8 Antibody, Biotin conjugated

Product Code	CSB-PA764872LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6RVD6
Immunogen	Recombinant Human Spermatogenesis-associated protein 8 protein (1-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenesis-associated protein 8 (Spermatogenesis-related protein 8), SPATA8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPATA8