





C6orf89 Antibody

Product Code	CSB-PA764939LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWU4
Immunogen	Recombinant Human Bombesin receptor-activated protein C6orf89 protein (120-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bombesin receptor-activated protein C6orf89 (Amfion), C6orf89, BRAP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C6orf89

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA764939LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA764939LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.