



C6orf89 Antibody, HRP conjugated

Product Code	CSB-PA764939LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWU4
Immunogen	Recombinant Human Bombesin receptor-activated protein C6orf89 protein (120-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bombesin receptor-activated protein C6orf89 (Amfion), C6orf89, BRAP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C6orf89
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.