



BTNL3 Antibody, Biotin conjugated

Product Code	CSB-PA764943LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXE8
Immunogen	Recombinant Human Butyrophilin-like protein 3 protein (18-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Butyrophilin-like protein 3 (Butyrophilin-like receptor), BTNL3, BTNLR COLF4100
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	BTNL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.