



SLC35D2 Antibody, HRP conjugated

Product Code	CSB-PA765728LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76EJ3
Immunogen	Recombinant Human UDP-N-acetylglucosamine/UDP-glucose/GDP-mannose transporter protein (75-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-N-acetylglucosamine/UDP-glucose/GDP-mannose transporter (Homolog of Fringe connection protein 1) (HFRC1) (SQV7-like protein) (SQV7L) (Solute carrier family 35 member D2) (UDP-galactose transporter-related protein 8) (UGTrel8), SLC35D2, HFRC UGTREL8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC35D2