



# SLC35D2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA765728LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q76EJ3
<b>Immunogen</b>	Recombinant Human UDP-N-acetylglucosamine/UDP-glucose/GDP-mannose transporter protein (75-146AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	UDP-N-acetylglucosamine/UDP-glucose/GDP-mannose transporter (Homolog of Fringe connection protein 1) (HFRC1) (SQV7-like protein) (SQV7L) (Solute carrier family 35 member D2) (UDP-galactose transporter-related protein 8) (UGTrel8), SLC35D2, HFRC UGTREL8
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC35D2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.