



CEP68 Antibody

Product Code	CSB-PA765741LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76N32
Immunogen	Recombinant Human Centrosomal protein of 68 kDa protein (280-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein of 68 kDa (Cep68), CEP68, KIAA0582
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEP68
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.