



TAS2R42 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA766403LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7RTR8 |
| Immunogen | Recombinant Human Taste receptor type 2 member 42 protein (149-187AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Taste receptor type 2 member 42 (T2R42) (Taste receptor type 2 member 55) (T2R55), TAS2R42, TAS2R55 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | TAS2R42 |