



BCDIN3D Antibody, HRP conjugated

Product Code	CSB-PA767206LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z5W3
Immunogen	Recombinant Human Pre-miRNA 5'-monophosphate methyltransferase protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-miRNA 5'-monophosphate methyltransferase (EC 2.1.1.-) (BCDIN3 domain-containing protein), BCDIN3D
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCDIN3D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.