





TADA2B Antibody

Product Code	CSB-PA768204LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TJ2
Immunogen	Recombinant Human Transcriptional adapter 2-beta protein (29-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional adapter 2-beta (ADA2-like protein beta) (ADA2-beta), TADA2B, ADA2B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TADA2B

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA768204LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA768204LA01HU at dilution of 1:100