



TSPAN33 Antibody, FITC conjugated

Product Code	CSB-PA768209LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UF1
Immunogen	Recombinant Human Tetraspanin-33 protein (118-235AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-33 (Tspan-33) (Penumbra) (hPen) (Proerythroblast new membrane), TSPAN33, PEN
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSPAN33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.