



# PAQR7 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA768224LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q86WK9
<b>Immunogen</b>	Recombinant Human Membrane progesterin receptor alpha protein (1-72AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Membrane progesterin receptor alpha (mPR alpha) (Membrane progesterone P4 receptor alpha) (Membrane progesterone receptor alpha) (Progesterone and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member VII), PAQR7, MRPA
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PAQR7
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.